SP-GRIPFLOW

Selective Cerebral Perfusion Catheter

Instructions for Use


SP-GRIPFLOW

Selective Cerebral Perfusion Catheter Instructions for use	EN	5
Cathéter de perfusion cérébrale sélective Mode d'emploi	FR	8
Selektiver zerebraler Perfusionskatheter Gebrauchsanweisung	DE	11
Catetere per perfusione cerebrale selettiva Istruzioni per l'uso	IT	14
Catéter de perfusión cerebral selectiva Instrucciones de uso	ES	17
Kateter för selektiv cerebral perfusion Bruksanvisningar	SV	20
Selectieve cerebrale perfusiekatheter Gebruiksaanwijzing	NL	23
Cateter de perfusão cerebral seletiva Instruções de utilização	PT	26
Selektiivinen aivoperfuusiokatetri Käyttöohjeita	FI	29
Kateter til selektiv cerebral perfusion Brugsanvisning	DA	32
Καθετήρας επιλεκτικής εγκεφαλικής αιμάτωσης Οδηγίες χρήσης	EL	35
Katétr pro selektivní mozkovou perfuzi Návod k použití	CS	38
Cewnik do selektywnej perfuzji mózgowej Instrukcja użycia	PL	41
Selektif Serebral Perfüzyon Kateteri Kullanım Talimatı	TR	44
Kateter for selektiv cerebral perfusjon Bruksinstruksjoner	NO	47
Selektivni kateter za cerebralnu perfuziju Upute za uporabu	HR	50

Device Description:

1. This product is a silicone catheter equipped with a silicone balloon to be used during surgery for aortic aneurysm.
2. During surgery for aortic arch aneurysm and dissecting aortic aneurysm, antegrade blood flow (perfusion) can be demonstrated for cerebroprotection through these catheters inserted into three branches of the aorta.
3. The balloon allows for an improved grip on/reduction of pulls from the target artery because of the rib on the balloon surface.
4. A catheter allows for adequate/optimizes perfusion flow rates.
5. Balloon compliance is referred to the Table 1.
6. There are two types of forms depending on the connector. Type I has a Luer Connector and Type II has a Taper connector.
7. The connector of each blood transmission line (MAIN) of this product can be connected to a blood circuit in a heart-lung machine equipped with a male luer-lock connector (Type I) or blood circuit tube (Type II) shown in the Table 2.

Figure 1:

A	Tip	G	Luer connector
B	Balloon	H	Cap
C	Marking	I	Cover tube (Red)
D	Shaft	J	Cover tube (Green)
E	Junction	K	Clamp (Red)
F	One-way valve	L	Taper connector

Table 1: Balloon Compliance

Injection Volume (mL)	Balloon Diameter (mm)		
	12F	14F	16F
0.0	6	7	8
1.0	12	11	11
2.0	16	16	15
3.0	18	18	18
3.5	19	19	19

*Maximum capacity of the balloon : 3.5 mL

Table 2: Connector Specification

Appropriate connector or tube	
Type I	A 6/100 tapered male luer-lock connector (ISO80369-7)
Type II	Blood circuit tube ID : 6.3mm (1/4 inch)

Indications for Use:

The SP-GRIPFLOW is indicated for selective cerebral perfusion of arteries originating from the aorta (brachiocephalic, left common carotid, and/or left subclavian).

Intended Purpose/Use:

The SP-GRIPFLOW is intended to be placed in the branches of the aortic arch for the cerebral perfusion during surgery of aorta.

Intended Users:

The SP-GRIPFLOW is intended to be used by physicians trained in interventional cardiovascular techniques and related surgical procedures. The devices are used in an operation room setting consisting of the appropriate equipment and personnel to support selective cerebral perfusion.

Intended Patient Population:

The SP-GRIPFLOW is intended for selective cerebral perfusion in subjects who will undergo an open surgical repair of the distal aortic arch.

Clinical benefit:

The SP-GRIPFLOW obstructs blood flow to the operative field (hemostatic techniques during surgery) and protects the brain (perfusion). In addition, the following benefits are anticipated for the device:

- The rib on the balloon surface allows for an improved grip and reduction of dislocation from the target artery.
- A wide catheter inner diameter allows for adequate perfusion flow rates.

Performance Characteristics:

- **Performance:** The SP-GRIPFLOW enables an adequate blood flow rate through the catheter and subsequent perfusion of the brain during open surgical repair of the aortic arch.
- **Safety:** The SP-GRIPFLOW is safe to use during open surgical repair of the distal aortic arch and does not cause an unacceptable rate of adverse events.

Contraindications:

1. This device is not to be used other than as intended use.
2. This device is not to be used for the patient with significant arteriosclerosis or calcified vascular walls.

Complications:

Possible complications include, but are not limited to, the following:
air embolism, cerebral infarction, cerebral ischemia, cerebral perfusion abnormalities, vessel damage

Warnings:

1. Do not reuse. This device is for single use only. Discard after one procedure. Structural integrity and / or function may be impaired through reuse or cleaning.
2. Do not re-sterilize. This device is disposable and sterilized by ethylene oxide gas.
3. Determine whether this product is appropriate for your patient by preoperative testing. [If the diameter of blood vessels is large, this product may not seal the blood vessels.]
4. Make sure to inflate the balloon with the sterile saline solution for balloon priming. [If air is injected into the balloon, the air may penetrate into the silicone balloon and cause an air embolism.]
5. The infusion amount into a balloon must be less than its maximum capacity based on clinical judgment. [Over inflation of the balloon can cause balloon rupture or blood vessel damage. (It is recommended that the surgeon check the balloon inflation status.)]
6. Determine the balloon placement site after undertaking diagnostic before and during surgery. [If there is an atheromatous plaque or mural thrombus at the placement site, it may become dislodged during placement, causing cerebral ischemia or cerebral infarction.]
7. Pay attention to the potential damage (perforation, dissection) to vessels caused by forcefully insertion of the catheter or over inflation of the balloon when placing this product.
8. When placing the catheter in arteries such as the brachiocephalic artery, place it in such a way that its tip does not over the vessel bifurcation. [If placing it over the vessel bifurcation, it may cause cerebral ischemia.]
9. If the balloon bursts, immediately stop perfusion (blood transmission), after no damage in the blood vessel is confirmed, exchange the catheter with a new one.

Precautions:

1. Do not use this catheter if the sterile package is damaged, wet or opened prior to use.
2. Do not use this catheter that has been damaged. [Damaged silicone products are significantly weakened.]
3. Do not grasp the catheter with forceps and the like. [If grasped with forceps and the like, the catheter may be damaged.]
4. Do not grasp the balloon with forceps and the like. [Grabbing the balloon with forceps and the like may damage it and cause it to burst.]
5. When inflating or deflating the balloon, make sure to insert the tip of the syringe all the way into the valve. [If the insertion is insufficient, the inside valve may not function, making it impossible to operate the balloon.]
6. Do not use the product if leakage and / or abnormal expansion is observed on the balloon during the inflation test.
7. During surgery, a spare (replacement) catheter is recommended to be prepared in case of unexpected catheter damage.
8. Make sure to connect the blood transmission line (MAIN) and pressure line (PRESSURE) correctly. [A green cover tube is attached to the distal part of the pressure line.]
9. Put the connector cap on the pressure measuring line (PRESSURE) properly after priming with the sterile heparinized saline solution if using the catheter without monitoring perfusion pressure based on clinical judgment. [Improper tightening of the connector cap may cause blood leakage.]
10. Connect each blood transmission line (MAIN) properly and check the connection status regularly.
11. Place (secure) the balloon at the appropriate placement site.
12. In case a vessel loop is used for securement, care should be taken to avoid physical damage to the catheter.
13. Insert the catheter and inflate the balloon while paying attention to the condition of the endothelium of the inserted blood vessel. (It is recommended that the surgeon check the balloon inflation status.)
14. Maintain an appropriate perfusion flow rate. (It is recommended to set the target flow rate at 10-13 mL/kg/min with brain perfusion temperature (Tympanic membrane temperature or nasal temperature) of 25°C) <Based on our clinical trial data>
15. Always confirm that there are no decreased volume of perfusion and/or abnormalities such as kinking, etc. during use.
16. Do not pull the shaft with a tension more than 15 N (1.5 kgf).
17. Remove the sterile saline solution in the balloon completely before performing intubation, extubation, and placement adjustment. [Incomplete deflation of the balloon may cause the burst and the damage to blood vessels.]
18. After use, dispose of this product and packaging based on the hospital, administrative and /or local government policy.
19. If user and/or patient notice any serious incident that has occurred in relation to the catheter, please report to the manufacturer and the competent authority.

Recommended Procedure:

1. Carefully remove the product from the sterile package, and check that the catheter is not damaged.
2. Prior to balloon inspection, inject and aspirate the sterile saline solution that is less than the maximum capacity for the purpose of removing air from the balloon and the balloon lumen (BALLOON). Repeat this procedure several times.
3. Inject the maximum volume of the sterile saline solution into the balloon, and check that there is no abnormality such as leakage or uneven inflation (balloon inspection)
4. Before inserting the catheter, the pressure measuring line (PRESSURE) and the blood transmission line (MAIN) must be adequately primed with the sterile heparinized saline solution.
5. After inserting the catheter into targeted blood vessels, inject the sterile saline solution into the balloon and place (secure) it. (Do not use air for balloon inflation)
6. If it is difficult to secure the catheter, it is possible to secure the catheter using a vessel loop.
7. It is recommended to connect the pressure measuring device to the pressure measuring line, and monitor the perfusion pressure at the catheter tip during perfusion.

Storage Conditions:

Keep the product dry and store in clean conditions, avoiding high temperature, humidity and direct sunlight.

Summary of Safety and Clinical Performance:

Information on the summary of safety and clinical performance (SSCP) is available at the following address:

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MedDevConsult GmbH
Fuji Systems Corporation, Shirakawa Plant